

Product Datasheet

Goat Rabbit IgG Secondary Antibody HRP Conjugated (orb43514)

Catalog Number	orb43514
Category	Antibodies
Description	Goat Anti-Rabbit IgG Secondary Antibody HRP Conjugated
Clonality	Polyclonal
Conjugation	HRP
Reactivity	Rabbit
Dilution range	IHC/ICC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:50,000 Western blot: 1:5,000 - 1:50,000
Application notes	IHC/ICC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:50,000 Western blot: 1:5,000 - 1:50,000
Antibody Type	Secondary Antibody
Storage	Store at -20°C/1 year
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)