

Product Datasheet

Human C-REACTIVE PROTEIN protein (orb435594)

Catalog Number	orb435594
Category	Proteins
Description	Native human C-REACTIVE PROTEIN protein
Target	C-REACTIVE PROTEIN
Reactivity	Human
Form/Appearance	Purified protein from human plasma - liquid
Concentration	1.0 mg/ml
Buffer/Preservatives	Purified protein from human plasma - liquid
Purity	>96% by SDS-PAGE
UniProt ID	P02741
Tested applications	ELISA
Specificity	C-REACTIVE PROTEIN
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the protein. Storage in frost-free freezers is not recommended.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558