

Product Datasheet

Human IgA protein (orb435626)

Catalog Number	orb435626
Category	Tools
Description	Native human IgA protein
Target	IgA
Reactivity	Human
Form/Appearance	Purified protein from human colostrum - liquid
Concentration	1.7 mg/ml
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃); Buffer: TRIS buffered saline
Purity	HPLC: >95%
UniProt ID	P01877, P01876
Tested applications	ELISA
Specificity	IgA
Storage	This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558