

Product Datasheet

Human RETINOL BINDING PROTEIN protein (orb435656)

Catalog Number	orb435656
Category	Proteins
Description	Native human RETINOL BINDING PROTEIN protein
Target	RETINOL BINDING PROTEIN
Reactivity	Human
Form/Appearance	Purified protein from the urine of patients with chronic renal tubular proteinuria - lyophilized
Concentration	1.0 mg/ml
Buffer/Preservatives	None present
Purity	>96% by SDS PAGE
UniProt ID	P02753
Tested applications	ELISA
Specificity	RETINOL BINDING PROTEIN
Storage	This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. After reconstitution store at -20°C. Storage in frost-free freezers is not recommended. Avoid repeated freezing and thawing as this may denature the protein.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)