

Product Datasheet

Human IL11 protein (orb435778)

Catalog Number	orb435778
Category	Proteins
Description	Recombinant of human IL11 protein
Target	IL-11
Reactivity	Human
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	100 ng/ul
Buffer/Preservatives	None present
Purity	>98% by SDS PAGE
UniProt ID	P20809
Tested applications	FA
Specificity	IL-11
Storage	Prior to reconstitution store at -20°C. Following reconstitution store at -20°C. Avoid repeated freezing and thawing as this may denature the protein. Should this product contain a precipitate we recommend microcentrifugation before use.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)