

## Product Datasheet

### Viral INFLUENZA B ANTIGEN protein (orb435903)

<b>Catalog Number</b>	orb435903
<b>Category</b>	Proteins
<b>Description</b>	Native viral INFLUENZA B ANTIGEN protein
<b>Target</b>	INFLUENZA B ANTIGEN
<b>Reactivity</b>	Virus
<b>Form/Appearance</b>	Inactivated Influenza B antigen - liquid
<b>Concentration</b>	Current, batch-specific concentration 0.91 mg/ml
<b>Buffer/Preservatives</b>	Preservative: None present; Buffer: Eagle's Minimum Essential Medium
<b>Purity</b>	Inactivated Pathogen
<b>Tested applications</b>	ELISA
<b>Specificity</b>	INFLUENZA B ANTIGEN
<b>Storage</b>	Store at -70°C. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the protein.
<b>Note</b>	For research use only
<b>Expiration Date</b>	6 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)