

Product Datasheet

GOAT ANTI MOUSE IgM:APC (HUMAN ADSORBED) (orb436164)

Catalog Number	orb436164
Category	Antibodies
Description	Goat polyclonal antibody against IgM conjugated to APC
Target	IgM
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	APC
Reactivity	Mouse
Form/Appearance	Purified Ig conjugated to Allophycocyanin (APC) - liquid
Concentration	0.5 mg/ml
Buffer/Preservatives	Purified Ig conjugated to Allophycocyanin (APC) - liquid
Immunogen	Mouse IgM paraproteins.
UniProt ID	P01872, P01873
Tested applications	ELISA, FACS, ICC, IF
Storage	+4°C DO NOT FREEZE
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)