

Product Datasheet

GOAT ANTI MOUSE IgM:APC (HUMAN ADSORBED) (orb436165)

Catalog Number	orb436165
Description	Goat polyclonal antibody against IgM conjugated to APC
Species/Host	Goat
Reactivity	Mouse
Conjugation	APC
Tested Applications	ELISA, FACS, ICC, IF
Immunogen	Mouse IgM paraproteins.
Target	IgM
Preservatives	Purified Ig conjugated to Allophycocyanin (APC) - liquid
Form/Appearance	Purified Ig conjugated to Allophycocyanin (APC) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P01873, P01872
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)