

Product Datasheet

GOAT ANTI RAT IgG (H/L):HRP (orb436200)

Catalog Number	orb436200
Category	Antibodies
Description	Goat polyclonal antibody against IgG conjugated to HRP
Target	IgG
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	HRP
Reactivity	Rat
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Buffer/Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Immunogen	Rat IgG, whole molecule.
UniProt ID	P20767, P20761, P20762, P20766, P20760, P01836, P01835, P20759
Tested applications	ELISA, IH-Fr, WB
Application notes	(H/L)
Storage	-20°C after reconstitution
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)