

Product Datasheet

HISTIDINE TAG Antibody (Biotin) (orb436318)

Catalog Number	orb436318
Category	Antibodies
Description	Goat polyclonal antibody against HISTIDINE TAG conjugated to Biotin
Target	HISTIDINE TAG
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	Biotin
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Buffer/Preservatives	Preservative: 0.09% Sodium Azide (NaN ₃) 0.2% Bovine Serum Albumin; Buffer: Phosphate buffered saline
Immunogen	Keyhole limpet haemocyanin (KLH) conjugated H-H-H-H-H-H
Tested applications	ELISA, WB
Specificity	HISTIDINE TAG
Antibody Type	Secondary Antibody
Storage	This product is shipped at ambient temperature. Store at +4°C. DO NOT FREEZE. This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)