

Product Datasheet

RAT ANTI MOUSE IgM HEAVY CHAIN (orb436347)

Catalog Number	orb436347
Category	Antibodies
Description	Rat monoclonal antibody to IgM
Target	IgM
Clonality	Monoclonal
Species/Host	Rat
Isotype	IgM
Conjugation	Unconjugated
Reactivity	Mouse
Form/Appearance	Purified IgM - liquid
Concentration	1 mg/ml
Buffer/Preservatives	Purified IgM - liquid
Purity	Purified
Immunogen	Polyclonal mouse IgM.
UniProt ID	P01872, P01873
Tested applications	ELISA
Application notes	HEAVY CHAIN
Clone Number	LO-MM-3

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)