

Product Datasheet

MOUSE ANTI RABBIT IgM (B CELL MARKER):RPE (orb436470)

Catalog Number	orb436470
Category	Antibodies
Description	Mouse monoclonal antibody against IgM conjugated to RPE
Target	IgM
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	RPE
Reactivity	Rabbit
Form/Appearance	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	Purified IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
UniProt ID	P03988, P04221
Tested applications	FC
Application notes	(B-CELL MARKER)
Clone Number	NRBM
Storage	+4°C DO NOT FREEZE
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)