

Product Datasheet

GOAT ANTI RABBIT IgG (Fc):HRP (orb436521)

Catalog Number	orb436521
Category	Antibodies
Description	Goat polyclonal antibody against IgG conjugated to HRP
Target	IgG
Clonality	Polyclonal
Species/Host	Goat
Isotype	Polyclonal IgG
Conjugation	HRP
Reactivity	Rabbit
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Buffer/Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
UniProt ID	P01870
Tested applications	ELISA, WB
Application notes	(Fc)
Storage	+4°C DO NOT FREEZE
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)