

Product Datasheet

RABBIT ANTI GOAT IgG (Fc):HRP (orb436524)

Description	Rabbit polyclonal antibody against IgG conjugated to HRP
Species/Host	Rabbit
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Target	IgG
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(Fc)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)