

Product Datasheet

GOAT ANTI RABBIT IgG (H/L):HRP (orb436525)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Rabbit IgG, whole molecule.
Target	IgG
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - lyophilized
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Uniprot ID	P01870, P01847, P01841, P01840, P01839, P01838
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)