

Product Datasheet

IgG Antibody (RPE) (orb436543)

Catalog Number	orb436543
Category	Antibodies
Description	Rabbit polyclonal antibody against IgG conjugated to RPE
Target	IgG
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Polyclonal IgG
Conjugation	RPE
Reactivity	Mouse
Form/Appearance	F(ab') ₂ fragment of IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Buffer/Preservatives	F(ab') ₂ fragment of IgG conjugated to R. Phycoerythrin (RPE) - lyophilized
Immunogen	Mouse IgG.
UniProt ID	P01863, P01869, P01868, P01864, P01867, P01865, P03987
Tested applications	FC
Specificity	IgG
Antibody Type	Secondary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

This product is shipped at ambient temperature. Prior to reconstitution store at +4°C. Following reconstitution store at +4°C. DO NOT FREEZE. This product should be stored undiluted. This product is photosensitive and should be protected from light. Should this product contain a precipitate we recommend microcentrifugation before use.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)