

Product Datasheet

IgG Antibody (orb436625)

Catalog Number	orb436625
Category	Antibodies
Description	Rabbit polyclonal antibody to IgG
Target	IgG
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	Polyclonal IgG
Conjugation	Unconjugated
Reactivity	Goat
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Buffer/Preservatives	Purified IgG - liquid
Purity	Purified
Immunogen	Normal goat IgG.
Tested applications	ELISA, ICC, IF, WB
Application notes	(H/L)
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)