

## Product Datasheet

### IgG IgA IgM Antibody (HRP) (orb436655)

|                             |                                                                                            |
|-----------------------------|--------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb436655                                                                                  |
| <b>Category</b>             | Antibodies                                                                                 |
| <b>Description</b>          | Goat polyclonal antibody against IgG IgA IgM conjugated to HRP                             |
| <b>Target</b>               | IgG IgA IgM                                                                                |
| <b>Clonality</b>            | Polyclonal                                                                                 |
| <b>Species/Host</b>         | Goat                                                                                       |
| <b>Isotype</b>              | Polyclonal IgG                                                                             |
| <b>Conjugation</b>          | HRP                                                                                        |
| <b>Reactivity</b>           | Mouse                                                                                      |
| <b>Form/Appearance</b>      | Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid                           |
| <b>Concentration</b>        | 1.0 mg/ml                                                                                  |
| <b>Buffer/Preservatives</b> | Preservative: 0.0095% MIT; Buffer: Phosphate buffered saline                               |
| <b>Immunogen</b>            | Mouse IgG purified from normal mouse serum, mouse IgA and IgM purified from ascitic fluid. |
| <b>UniProt ID</b>           | <b>P03987, P01878, P01869, P01864, P01865, P01868, P01867, P01873, P01863, P01872</b>      |
| <b>Tested applications</b>  | ELISA                                                                                      |
| <b>Specificity</b>          | IgG IgA IgM                                                                                |
| <b>Antibody Type</b>        | Secondary Antibody                                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Storage**

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C. Avoid repeated freezing and thawing as this may denature the antibody. Storage in frost-free freezers is not recommended.

**Note**

For research use only

**Expiration Date**

12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)