

Product Datasheet

Human IL1e ELISA Kit (orb437240)

Catalog Number	orb437240
Category	Assays and Kits
Description	Human Interleukin 1 Epsilon (IL1e) Sandwich ELISA Kit
Reactivity	Human
Form/Appearance	Traditional
UniProt ID	Q9UHA7
Assay Type	Sandwich
Assay Time	4.5 hours
Range	15.625-1000pg/mL
Sensitivity	5.1pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Entrez	27179

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)