

Product Datasheet

Chicken CCKAR ELISA Kit (orb438438)

Catalog Number	orb438438
Description	Chicken Cholecystokinin A Receptor (CCKAR) Sandwich ELISA Kit
Reactivity	Gallus
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: tissue homogenates and other biological fluids; sandwich: Experimental Method
Sample Types	tissue homogenates and other biological fluids.
Assay Time	4.5 hours
Uniprot ID	E5RVJ7
Sensitivity	0.179ng/mL
Entrez	422801
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)