

Product Datasheet

Anti-GAPDHS Purified (orb43877)


Description	Mouse Monoclonal to Intra-Acrosomal Protein.
Reactivity	Human, Porcine
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	Freshly ejaculated human sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid and subsequently used for immunization.
Target	GAPDHS
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Immunocytochemistry: Recommended dilution: 10 µg/ml; membrane permeabilization (acetone) is essential. Flow cytometry: Recommended dilution: 3-12 µg/ml. Intraacrosomal staining.
Isotype	Mouse IgM
Clonality	Monoclonal
Clone Number	Hs-8
Antibody Type	Primary Antibody
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Uniprot ID **O14556****Entrez** **26330****Expiration Date** 12 months from date of receipt.

Separation of human GAPDHS positive cells (red-filled) from GAPDHS negative cells (black-dashed) in flow cytometry analysis (intracellular staining) of human sperm cells stained using anti-GAPDHS (Hs-8) purified antibody (concentration in sample 7.5 µg/ml) GAM APC.



Flow cytometry intracellular staining pattern of human sperm cells stained using anti-GAPDHS (Hs-8) purified antibody (concentration in sample 7.5 µg/ml) GAM APC.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)