

Product Datasheet

General PGF2a ELISA Kit (orb441534)

Catalog Number	orb441534
Category	Assays and Kits
Description	General Prostaglandin F2 Alpha (PGF2a) competitive inhibition ELISA Kit
Reactivity	All
Form/Appearance	Traditional
Assay Type	Competitive
Assay Time	4.5 hours
Range	31.25-2000pg/mL
Sensitivity	7.85pg/mL
Sample Types	serum, plasma and other biological fluids.
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)