

## Product Datasheet

### Human UBE2C ELISA Kit (orb442800)

<b>Catalog Number</b>	orb442800
<b>Description</b>	Human Ubiquitin Conjugating Enzyme E2C (UBE2C) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	Sample Type: tissue homogenates, cell lysates and other biological fluids; sandwich: Experimental Method
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Assay Time</b>	4.5 hours
<b>Uniprot ID</b>	<b>Q5TZN3</b>
<b>Sensitivity</b>	0.053ng/mL
<b>Entrez</b>	<b>11065</b>
<b>Expiration Date</b>	Please enquire.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)