

Product Datasheet

Mouse UCN3 ELISA Kit (orb442821)

Catalog Number	orb442821
Description	Mouse Urocortin 3 (UCN3) Sandwich ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	Sample Type: serum, plasma and other biological fluids; sandwich: Experimental Method
Sample Types	serum, plasma and other biological fluids.
Assay Time	4.5 hours
Uniprot ID	Q924A4
Sensitivity	4.3pg/mL
Entrez	83428
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)