

Product Datasheet

Mouse CRLF2/TSLP R ELISA Kit (orb443256)

| | |
|----------------------------|---|
| Catalog Number | orb443256 |
| Category | Assays and Kits |
| Description | Mouse CRLF2/TSLP R ELISA Kit |
| Target | Cytokine receptor-like factor 2 |
| Reactivity | Mouse |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: A20-L233 |
| UniProt ID | Q8CII9 |
| Tested applications | ELISA |
| Specificity | Natural and recombinant mouse CRLF2 |
| Cross Reactivity | There is no detectable cross-reactivity with other relevant proteins. |
| Range | 125 pg/ml - 8,000 pg/ml |
| Sensitivity | 10 pg/ml |
| Sample Types | cell culture supernatants, serum and plasma (heparin, EDTA). |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Expiration Date | 6 months from date of receipt. |

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Mouse CRLF2/TSLP R PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)