

Product Datasheet

Anti-Cytokeratin 10 Purified (orb44573)

Description	Mouse Monoclonal to Cytokeratin 10.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, IP
Immunogen	Cytoskeleton preparation extracted from human epidermis by detergent/high salt extraction.
Target	Cytokeratin 10
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: skin, pretreatment: trypsin/temperature 37 °C.
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	VIK-10
Antibody Type	Primary Antibody
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P13645

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Entrez**3858****Expiration Date**

12 months from date of receipt.



Immunohistochemistry staining of skin basaliom (paraffin-embedded sections) with anti-human Cytokeratin 10 (VIK-10).

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44%201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)