

## Product Datasheet

### MCL1 Antibody (FITC) (orb46201)

|                             |                                                                                                                                                                                                                                                      |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb46201                                                                                                                                                                                                                                             |
| <b>Category</b>             | Antibodies                                                                                                                                                                                                                                           |
| <b>Description</b>          | This MCL1 Antibody, FITC conjugated is an FITC polyclonal product. It tagets MCL1 using a Recombinant Human Induced myeloid leukemia cell differentiation protein Mcl-1 protein (19-244AA) as the immunogen. Purification: >95%, Protein G purified. |
| <b>Target</b>               | MCL1                                                                                                                                                                                                                                                 |
| <b>Clonality</b>            | Polyclonal                                                                                                                                                                                                                                           |
| <b>Species/Host</b>         | Rabbit                                                                                                                                                                                                                                               |
| <b>Isotype</b>              | IgG                                                                                                                                                                                                                                                  |
| <b>Conjugation</b>          | FITC                                                                                                                                                                                                                                                 |
| <b>Reactivity</b>           | Human                                                                                                                                                                                                                                                |
| <b>Form/Appearance</b>      | Liquid                                                                                                                                                                                                                                               |
| <b>Buffer/Preservatives</b> | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                                                                                                                                                                                          |
| <b>Purification</b>         | >95%, Protein G purified                                                                                                                                                                                                                             |
| <b>Immunogen</b>            | Recombinant Human Induced myeloid leukemia cell differentiation protein Mcl-1 protein (19-244AA)                                                                                                                                                     |
| <b>UniProt ID</b>           | <b>Q07820</b>                                                                                                                                                                                                                                        |
| <b>Dilution range</b>       | ELISA                                                                                                                                                                                                                                                |
| <b>Storage</b>              | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                        |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Note**

For research use only

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)