

## Product Datasheet

# Phospho-PFKFB3/PFK2 (Ser467) Rabbit Polyclonal Antibody (HRP) (orb480027)

|                             |  |
|-----------------------------|--|
| <b>Catalog Number</b>       | orb480027  |
| <b>Category</b>             | Antibodies   |
| <b>Description</b>          | Phospho-PFKFB3/PFK2 (Ser467) Rabbit Polyclonal Antibody (HRP) is a HRP conjugated antibody targeting PFKFB3. This antibody is suitable for IHC-Fr, IHC-P, WB. It exhibits reactivity with Human, Mouse, Rat samples. |
| <b>Target</b>               | PFKFB3   |
| <b>Clonality</b>            | Polyclonal   |
| <b>Species/Host</b>         | Rabbit   |
| <b>Isotype</b>              | IgG  |
| <b>Conjugation</b>          | HRP  |
| <b>Reactivity</b>           | Human, Mouse, Rat  |
| <b>Predicted Reactivity</b> | Bovine, Equine, Gallus, Porcine, Rabbit  |
| <b>Concentration</b>        | 1mg/ml   |
| <b>Buffer/Preservatives</b> | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.   |
| <b>Purification</b>         | Affinity purified by Protein A   |
| <b>Immunogen</b>            | KLH conjugated synthesised phosphopeptide derived from human PFK2 around the phosphorylation site of Ser467 LA(p-S)PE  |
| <b>Tested applications</b>  | IHC-Fr, IHC-P, WB  |

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

|                        |   |
|------------------------|---|
| <b>Dilution range</b>  | WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500                                 |
| <b>Antibody Type</b>   | Primary Antibody  |
| <b>Storage</b>         | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |

---

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)