

Product Datasheet

GFP Mouse mAb, PE conjugated (orb490698)

Catalog Number	orb490698
Category	Antibodies
Description	Mouse monoclonal antibody against GFP-Tag conjugated to PE
Target	GFP
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG
Conjugation	PE
Reactivity	All
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Buffer/Preservatives	Lyophilized or Liquid
Purification	Affinity purified by Protein G
Immunogen	Recombinant Green Fluorescence Protein
Tested applications	FC, ICC, IF
Dilution range	IF=1:200-1000, Flow-Cyt=1ug/Test, ICC/IF=1:100-500
Antibody Type	Primary Antibody
Clone Number	B6I8
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)