

Product Datasheet

Human SUMF1 protein (orb49670)

Catalog Number	orb49670
Category	Proteins
Description	Recombinant Human Sulfatase Modifying Factor 1 is produced by our Mammalian expression system and the target gene encoding Ser34-Asp374 is expressed with a 6His tag at the C-terminus.
Tag	C-His
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, 2mM CaCl ₂ , 10% Glycerol, pH 7.5.
UniProt ID	Q8NBK3
Application notes	6His tag at the C-terminus.
Expression System	Human Cells
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Biological Origin	Human
Expression Region	Ser34-Asp374
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Dry Ice Shipping	Please note: This product requires shipment on dry ice. A dry ice surcharge will apply.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558