

Product Datasheet

UCHL1+3 Rabbit pAb, PE conjugated (orb499563)

Catalog Number	orb499563
Category	Antibodies
Description	UCHL1+3 Rabbit pAb, PE conjugated is a PE conjugated antibody targeting UCHL1. This antibody is suitable for ICC. It exhibits reactivity with Human, Mouse, Rat samples.
Target	UCHL1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PE
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Porcine, Rabbit, Sheep
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Buffer/Preservatives	Lyophilized or Liquid
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated synthetic peptide derived from human UCHL1+3 (141-223/223aa)
Tested applications	ICC
Dilution range	ICC/IF=1:50-200

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)