

Product Datasheet

Alpha Tubulin (Acetyl Lys40) Mouse Monoclonal Antibody (orb499956)

Description	Alpha Tubulin (Acetyl Lys40) Mouse Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised acetylpeptide derived from human alpha Tubulin around the acetylation site of Acetyl Lys40 SD(Ac-K)TI
Target	TUBA4A
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Clone Number	6B5
Antibody Type	Primary Antibody
MW	50 kDa

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Uniprot ID


P68366

Dilution Range


WB=1:500-2000

Expiration Date

12 months from date of receipt.



Sample: HepG2 (Human) Cell Lysate at 30 ug, HeLa (Human) Cell Lysate at 30 ug, NIH/3T3 (Mouse) Cell Lysate at 30 ug, A431 (Human) Cell Lysate at 30 ug, Primary: Anti-alpha Tubulin (Acetyl Lys40) (orb499956) at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution, Predicted band size: 52 kD, Observed band size: 53 kD.



Sample: Lane 1: A549 (Human) Cell Lysate at 30 ug, Lane 2: LOVO (Human) Cell Lysate at 30 ug, Lane 3: A431 (Human) Cell Lysate at 30 ug, Lane 4: HeLa (Human) Cell Lysate at 30 ug, Primary: Anti-alpha Tubulin (Acetyl Lys40) (orb499956) at 1/1000 dilution, Secondary: IRDye800CW Goat Anti-Mouse IgG at 1/20000 dilution, Predicted band size: 52 kD, Observed band size: 52 kD.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)