

Product Datasheet

Human MCP-1 / CCL2 ELISA Kit (orb50069)

Catalog Number	orb50069
Category	Assays and Kits
Description	Human MCP-1 / CCL2 ELISA Kit
Target	C-C motif chemokine 2
Reactivity	Human
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-T99
UniProt ID	P13500
Tested applications	ELISA
Specificity	Natural and recombinant human MCP-1
Cross Reactivity	There is cross-reactivity with human Eotaxin, MCP-31%.
Range	15.6 pg/ml - 1,000 pg/ml
Sensitivity	1 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human MCP-1 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)