

Product Datasheet

Human sTNFsR II ELISA Kit (orb50114)

Catalog Number	orb50114
Category	Assays and Kits
Description	Human sTNFsR II ELISA Kit
Target	Tumor necrosis factor receptor superfamily member 1B
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: L23-D257
UniProt ID	P20333
MW	48291 Da
Tested applications	ELISA
Specificity	Natural and recombinant human sTNFsRII
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	7.8 pg/ml - 500 pg/ml
Sensitivity	2 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Human sTNFsR II ELISA Kit

Human sTNFsR II PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)