

Product Datasheet

Human IL-18/IL-1F4/Interleukin-18 ELISA Kit (orb50153)

Catalog Number	orb50153
Category	Assays and Kits
Description	Human IL-18/IL-1F4/Interleukin-18 ELISA Kit
Target	Interleukin-18
Reactivity	Human
UniProt ID	Q14116
MW	22326 Da
Tested applications	ELISA
Specificity	Natural and recombinant human IL-18
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Standard	Expression system for standard: E.coli; Immunogen sequence: Y33-D197
Range	31.2 pg/ml - 2,000 pg/ml
Sensitivity	1 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Storage	Store at 4°C for 6 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 6 months from date of receipt.

Human IL-18 ELISA Kit



Human IL-18/IL-1F4/Interleukin-18 ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)