

Product Datasheet

phospho-IGF1R (Tyr1165 + Tyr1166) Rabbit pAb, HRP conjugated (orb503127)

Catalog Number	orb503127
Category	Antibodies
Description	Rabbit polyclonal antibody against IGF1R (phospho-Tyr1165/Tyr1166) conjugated to HRP
Target	IGF1R
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	HRP
Reactivity	Mouse, Rat
Predicted Reactivity	Bovine, Canine, Gallus, Human, Porcine
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Buffer/Preservatives	Lyophilized or Liquid
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IGF1R around the phosphorylation site of Tyr1165+Tyr1166 TD(p-Y)(p-Y)RK
UniProt ID	P08069

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

MW	150kDa
Tested applications	ELISA
Dilution range	ELISA=1:5000-1000
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)