

Product Datasheet

HRP Goat Mouse IgG Fc Antibody (orb511055)

Catalog Number	orb511055
Category	Antibodies
Description	Goat Anti-Mouse IgG Fc recognizes mouse IgG F(c) fragment. This HRP-conjugated antibody was purified using antigen affinity techniques and is suitable for ELISA and Western blotting procedures. The IgG component of this antibody conjugate is specific for normal mouse IgG and is non-reactive with IgM or IgE serum proteins. The IgG-HRP conjugate contains between 2 and 4 HRP molecules per IgG molecule.
Conjugation	Unconjugated
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)