

Product Datasheet

TRIM26 Antibody, FITC conjugated (orb516646)

Catalog Number	orb516646
Category	Antibodies
Description	This TRIM26 Antibody, FITC conjugated is an FITC polyclonal product. It targets TRIM26 using a Recombinant Human Tripartite motif-containing protein 26 protein (151-270AA) as the immunogen. Purification: >95%, Protein G purified.
Target	TRIM26
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	FITC
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Tripartite motif-containing protein 26 protein (151-270AA)
UniProt ID	Q12899
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)