

Product Datasheet

Mouse SerpinE2 ELISA Kit (orb526840)

Catalog Number	orb526840
Category	Assays and Kits
Description	Mouse SerpinE2 ELISA Kit
Target	Glia-derived nexin
Reactivity	Mouse
Immunogen	Expression system for standard: NS0; Immunogen sequence: S20-P397
UniProt ID	Q07235
Tested applications	ELISA
Specificity	Natural and recombinant mouse SerpinE2
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	0.78 ng/ml - 50 ng/ml
Sensitivity	0.1 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate) ,urine.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Mouse SerpinE2 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)