

Product Datasheet

PBS Buffer pH 7.6 - 10x conc. (orb532232)

Catalog Number	orb532232
Description	A 1x concentration of PBS (Phosphate buffered saline). Buffer contains: 8.01 g/l NaCl 0.2 g/l KCl 1.78 g/l Na 2HPO 4·2H2O 0.27 g/l KH 2PO4. The Buffer is used for cell culture applications since it is isotonic and non-toxic to cells.
Form/Appearance	liquid; pH: 7.6
Concentration	10x
Storage	store at 4 °C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)