

Product Datasheet

Tris-HCl Buffer 1 M pH 8.0 (orb532242)

Description	Tris-HCl Buffer is included in PCR reaction buffer, used for in vitro cell culture, enzyme assays and electrophoretic applications.
Form/Appearance	liquid; pH: 8.0 ± 0.1 (22 °C)
Concentration	1 M
Storage	store at 4 °C
Note	For research use only
Formula	$C_4H_{11}NO_3 \cdot HCl$
MW	Theoretical MW: 157.60 g/mol
CAS Number	1185-53-1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)