

Product Datasheet

5-Propargylamino-dCTP-Cy5 (orb533141)

Catalog Number	orb533141
Category	Tools
Description	5-Propargylamino-dCTP-Cy5 is recommended for direct enzymatic labeling of DNA/cDNA e.g. by PCR and Nick Translation. It is incorporated as substitute for its natural counterpart dCTP. The resulting Dye-labeled DNA/cDNA probes are ideally suited for fluorescence hybridization applications such as FISH or microarray-based gene expression profiling. Optimal substrate properties and thus labeling efficiency is ensured by an optimized linker attached to the C5 position of cytidine. Recommended Propargylamino-dCTP-Cy5/dCTP ratio for PCR and Nick Translation: 30-50% Propargylamino-dCTP-Cy5/70-50% dCTP.
Form/Appearance	solution in 10 mM Tris-HCl; Color: blue; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1159.01 g/mol (free acid); Detected MW: 1158.23 g/mol (free acid)
Application notes	Incorporation into DNA/cDNA by: PCR with Taq polymerase in-house data; Nick Translation with DNase I/DNA Polymerase in-house data. Spectroscopic Propertie: λ _{exc} 649 nm, λ _{em} 670 nm, ε 250.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₄₅ H ₅₇ N ₆ O ₂₀ P ₃ S ₂
SMILES	<chem>N(CC#Cc1cn(c(=O)nc1N)[C@@H]1O[C@H](COP(=O)(O)OP(=O)(O)OP(=O)(O)O)[C@H](C1O)C(=O)CCCCCN1c2ccc(S(=O)(=O)O)cc2C(C)(C)/C1=C\C=C\C=C\C1=*(c2ccc(cc2C1(C)C)S(=O)(=O)[O-])CC</chem>
Storage	store at -20 °C

Note For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558