

Product Datasheet

Crystal Taq Master (2x) (orb533329)

Catalog Number	orb533329
Description	Crystal PCR Master is a 2 x conc. ready-to-use master mix recommended for routine PCR applications (up to 4 kb fragment length), high throughput PCR or genotyping. It contains all reagents required for PCR (except template and primer) in a well-balanced ratio to ensure high specificity and minimal by-product formation in almost all PCR applications without the need of additional optimization steps. The mix guarantees robust and reliable amplification results with a minimum of pipetting steps, saves time and reduces the risk of contaminations. The total PCR assay volume is freely adaptable to individual protocols or the requirements of automated pipetting systems.
Form/Appearance	liquid
Concentration	2x conc
Storage	store at -20 °C
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)