

## Product Datasheet

### Blood DNA Preparation - Solution Kit (orb533439)

<b>Catalog Number</b>	orb533439
<b>Category</b>	Assays and Kits
<b>Description</b>	Blood DNA Preparation Kit is designed for convenient and fast isolation of genomic DNA from whole blood samples. The solution based system minimizes DNA fragmentation that may be problematic in spin-column/filtration based methods. Because phenol or chloroform is not used it is safe and does not produce any harmful waste. Solution based genomic DNA purification kits guarantee minimal DNA fragmentation and yield DNA sized up to 150 kb.
<b>Storage</b>	store at ambient temperature
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

---

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713-2847  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)