

Product Datasheet

Hemoglobin Colorimetric Microplate Assay Kit (orb545608)

Catalog Number	orb545608
Category	Assays and Kits
Description	Hemoglobin Colorimetric Microplate Assay Kit which is used to detect and quantify the hemoglobin content in serum, with the detection range is 1 µg/ml - 100 µg/ml.
Reactivity	All
Tested applications	FA
Application notes	When testing hemoglobin in red blood cell samples, it is recommended to first lyse the red blood cells, and then collect the supernatant through centrifugation for testing.
Assay Type	Detection
Range	1 µg/ml - 100 µg/ml
Sample Types	Serum
CAS Number	9008-02-0
Storage	Shipped and store at 4°C for 6 months.
Hazard Information	Non-Toxic
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Representative Standard Curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)