

Product Datasheet

Human BMP-2 Quick ELISA Kit (orb546890)

Catalog Number	orb546890
Category	Assays and Kits
Description	Human BMP-2 Quick ELISA Kit
Target	Bone morphogenetic protein 2
Reactivity	Human
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q283-R396
UniProt ID	C8C060
MW	44702 Da
Tested applications	ELISA
Specificity	Natural and recombinant human BMP-2
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	31.2 pg/ml - 2,000 pg/ml
Sensitivity	2 pg/ml
Sample Types	bone tissue, cell culture supernatants and serum.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Human BMP-2 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)