

Product Datasheet

Human IL-5/Interleukin-5 Quick ELISA Kit (orb546902)

Catalog Number	orb546902
Category	Assays and Kits
Description	Human IL-5/Interleukin-5 Quick ELISA Kit
Target	Interleukin-5
Reactivity	Human
Immunogen	Expression system for standard: NS0; Immunogen sequence: I20-S134
UniProt ID	P05113
MW	15238 Da
Tested applications	ELISA
Specificity	Natural and recombinant human IL-5
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Range	7.8 pg/ml - 500 pg/ml
Sensitivity	2 pg/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Expiration Date 6 months from date of receipt.



Human IL-5 PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)