

Product Datasheet

Mouse OPG (TNFRSF11B) / Osteoprotegerin Quick ELISA Kit (orb546906)

Catalog Number	orb546906
Category	Assays and Kits
Description	Mouse OPG (TNFRSF11B) / Osteoprotegerin Quick ELISA Kit
Target	Tumor necrosis factor receptor superfamily member 11B
Reactivity	Mouse
Immunogen	Expression system for standard: NS0; Immunogen sequence: E22-L401
UniProt ID	O08712
MW	45923 Da
Tested applications	ELISA
Specificity	Natural and recombinant mouse OPG
Cross Reactivity	There is no detectable cross-reactivity with other relevant proteins.
Standard	10 ng/tube
Range	93.7 pg/ml - 6,000 pg/ml
Sensitivity	5 pg/ml
Sample Types	cell culture supernatants, serum and plasma (EDTA).
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only



Mouse OPG (TNFRSF11B) PicoKine ELISA Kit standard curve

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.comPhone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)**Biorbyt LLC**

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.comPhone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)