

Product Datasheet

GRP78 BiP HSPA5 Rabbit Monoclonal Antibody (orb547543)

Catalog Number	orb547543
Category	Antibodies
Description	GRP78 BiP HSPA5 Rabbit Monoclonal Antibody
Clonality	Monoclonal
Species/Host	Rabbit
Isotype	Rabbit IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Liquid
Concentration	0.5mg/ml
Purification	Affinity-chromatography
Immunogen	A synthesized peptide derived from human GRP78 BiP
UniProt ID	P11021
MW	28 kDa
Tested applications	ICC, IF, IHC, WB
Application notes	WB 1:500-1:2000IHC 1:50-1:200ICC/IF 1:50-1:200. Restore with deionized water (or equivalent) for reconstitution volume of 1.0 mL

Antibody Type	Primary Antibody
----------------------	------------------

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Clone Number BOO-8

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


Note For research use only



Immunohistochemical analysis of paraffin-embedded human pancreas, using GRP78 BiP Antibody.



Immunofluorescent analysis of HepG2 cells, using GRP78 BiP Antibody.



WB analysis of GRP78 BiP expression in (1) LnCaP Cell;(2) HepG2 Cell.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558