

Product Datasheet

Guinea pig IL4 ELISA Kit (Ready to Use) (orb548746)

Catalog Number	orb548746
Category	Assays and Kits
Description	Guinea pig Interleukin 4 (IL4) Sandwich ELISA KitReady to Use: Pre-diluted detection reagents for convenient and simple assays.
Reactivity	Guinea pig
Form/Appearance	Ready to Use
UniProt ID	G0WM66
Assay Type	Sandwich
Assay Time	3 hours
Range	1.56-100pg/mL
Sensitivity	0.57pg/mL
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids
Storage	Store at -20°C
Note	For research use only
Entrez	100720403

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)